



**HCAI / AMR National Team**  
**Technical Document**

**Environmental Testing for Carbapenemase Producing Enterobacterales (CPE)**

Document Type	Technical Document	Document developed by	HSE HCAI/AMR Clinical Programme
Approval Date	7 June 2018	Document author	HCAI/AMR Team
Document reference number	HCAI/AMR T001	Document approved by	Prof Martin Cormican MCRN 011105
Version number	1.0	Responsibility for implementation	Guide technical document if required by IPCT/Microbiology
Revision Date	12 months	Responsibility for review and audit	Prof Martin Cormican MCRN 011105
Draft or Final document	Final document	Pages	4

**Issued by M Cormican (MCRN 011105) National Lead HCAI & AMR July 2018**



**Antimicrobial Resistance  
and Infection Control Team**

## Introduction

Enterobacteriales can persist in suitable environmental niches for extended periods. Environmental contamination has emerged as potential reservoir of CPE that may contribute to colonisation or infection of patients. In some cases Infection Prevention and Control (IPC) teams may consider it appropriate to test clinical environments for CPE. Sensitivity of detection of CPE from the environment is dependent on the sampling method and the analytical method. This document describes an approach to testing for environmental contamination with CPE that has been applied in one hospital with a significant yield of positive results. The approach is outlined here for consideration for use by other IPC teams and laboratories. The method has not been comprehensively compared with other approaches some of which may be equally or more effective.

## Scope

This document is intended for use by IPC teams and laboratory practitioners.

## Approach to Testing

Aseptic technique is essential for the following procedure.

### Swab method for environmental sampling<sup>1</sup>

**Recommended:** Flocked nylon swabs or macrofoam swabs for maximum recovery of organisms<sup>2</sup>.

1. Perform alcohol rub hand hygiene, put on gloves.
2. Remove a swab from the sterile wrapping and moisten the tip by immersing it in a tube containing dilution liquid or neutralizing fluid.
3. Press the tip of the swab against the wall of the tube to remove excess liquid.
4. Place the sterile template on the surface to be investigated, place the tip of the swab on the surface inside the template and streak the area within, whilst rotating the swab between the thumb and forefinger in two directions at right angles to each other.
5. Put the swab into a tube containing 10 ml of dilution liquid (buffered saline/nutrient broth/ buffered peptone saline/ Brain Heart Infusion Broth) and aseptically break or cut off the stick.
6. After sample collection remove gloves and perform hand hygiene. Label the sample and submit for laboratory processing (*skip to step 7 "Laboratory Processing"*)

### Sponge method for environmental sampling<sup>1</sup>

1. Perform alcohol rub hand hygiene, put on gloves.
2. Open the plastic bag or container containing the sponge.
3. Remove the sponge aseptically from its container with a sterile forceps or a sterile gloved hand. Moisten the cloth or sponge with sufficient quantity of diluent (without excess).

**Note:** Commercially available sponges are pre moistened with neutralizing buffer.

4. Sample the chosen surface in two perpendicular directions, changing the face of the cloth or sponge. Return the sponge to the sterile container and close.
5. After sample collection remove gloves and perform hand hygiene. Label the sample and submit for laboratory processing.

## Laboratory Processing


6. Add 100mls of diluents (buffered saline/nutrient broth/buffered peptone saline) to the transport container containing the sponge. Replace the lid and mix thoroughly by shaking/vortex for 30 seconds.
7. Incubate the container at 37°C in aerobic conditions for 24 hours.
8. Re-suspend broth (Vortex in safety cabinet) and inoculate the Chromogenic CPE agar with 10µls of broth (avoiding aerosols), streak out for discreet colonies.

**Recommended:** Use a selective Chromogenic agar that is optimised for types of CPE that are common in Ireland including OXA-48.

9. Incubate the plate at 37°C in aerobic conditions for 24 hours.
10. Read plate for suspect colonies and identify all isolates e.g. MALDI.
11. Test isolates of Enterobacterales for CPE enzymes/CPE encoding genes by molecular analysis/ immuno-chromatography assay.
12. Send confirmed or suspect CPE to the National CPE Reference Laboratory Service.

**Note:** Bacteria of orders other than Enterobacterales may be recovered including *Acinetobacter spp.* and *Pseudomonas spp.* The relevance of these non- Enterobacterales will depend on the infection control issue under investigation but is not addressed further in this note on methods.

## Appendix 1: Flow Chart of Environmental Testing for Carbapenemase Producing Enterobacterales



## References:

1. BS ISO 18593:2004 Microbiology of food and animal feeding stuffs- Horizontal methods for sampling techniques from surfaces using contact plates and swabs.
2. S. Galvin, A. Dolan, O. Cahill, S. Daniels, H. Humphreys. Microbial monitoring of the hospital environment: why and how? *Journal of Hospital Infection* 2012; **82**: 143-151.
3. Acknowledgement: This document is based on work performed by B. Hananahoe, Surveillance Scientist, Galway University Hospital.